

## MiroValve® AUTOM. CATHETER VALVE FOR UROLOGY

### MiroValve® technical data

product code 089805


box 10 pcs

- ✓ can't forget to open or close the valve
- ✓ automatically open once connected to the bag
- ✓ the valve closes automatically when bag is removed
- ✓ easy operation with one hand
- ✓ no dripping of urine
- ✓ applicable with all common catheters and bags
- ✓ permanent automatic indwelling by connection to the bag
- ✓ good for bladder training
- ✓ easy spray disinfection
- ✓ coloured operating elements
- ✓ valve system will be changed together with catheter
- ✓ can be used for biliary drainage too

All products and packaging are latex and pvc free.  
A product video is provided here:  
<http://www.youtube.com/watch?v=A88emF907pE>



## MiroValve®



MiroValve® Automatic Catheter Valve